


OEM X-Ray Tubes for XRF

Features:

- **4-50 kV**
- **2-10 W**
- **Anode materials: Ag, Au, Rh, W, and others**
- **Compact**
- **Self-shielded**
- **End-window X-ray tube design**
- **Custom packaging for OEM applications**
- **Ideal for hand-held or benchtop XRF instruments**



OEM X-ray tubes: Custom packaging to meet the unique demands of the OEM.



Complete OEM systems for XRF: Includes detector, preamplifier, digital processor, power supplies, X-ray tube, and software (acquisition, control, interface).

NOTE: The OEM X-ray tubes are not the same as the Amptek Mini-X. The Mini-X has a USB interface to control the voltage and current through PC software. The OEM versions are controlled by user supplied analog voltages.

X-Ray Tube Contact

Newton Scientific, Inc.
255 Bent Street, Cambridge, MA 02141, USA
TEL: +1 (617) 354-9469, FAX: +1 (617) 354-9479
info@newtonscientificinc.com
www.newtonscientificinc.com

X-Ray Detector Contact

Amptek, Inc.
14 DeAngelo Drive, Bedford, MA. 01730, USA
TEL: +1 (781) 275-2242, FAX: +1 (781) 275-3470
sales@amptek.com
www.amptek.com

NSI

